



Resting-State fMRI: ICA and Dual Regression

FSL Course 2024

19 June, Osaka, Japan



Resting state fMRI and ICA


- Introduction to resting state
- Independent Component Analysis
- Single-subject ICA
- Multi-subject ICA
- Dual regression



Why resting state fMRI?

Energy consumption in the brain

- Brain < 2% body weight but consumes ~20% of total energy
- Estimated 60-80% of this energy used to support communication between cells
- fMRI provides a window to brain activity
- Task-evoked activity accounts for ~1%




Finn 2021


Resting state -> intrinsic functional brain organisation

Why study the brain at rest?

- Understand the inherent functional organisation of the brain
- Clinical/ cognitive biomarker
- Pragmatic benefits: can be done in any population, with relatively little setup and expertise required
- Localisation vs connectivity



Biswal et al (1995)



Sheline et al (2010)



Overview of Resting State Analysis



Principles of resting state analysis


- Many different methods available for analysis
- All have one assumption in common:
 - The definition of functional connectivity is based on a statistical dependency between timeseries
- Differences between methods lie in the way these similarities are estimated and/or represented

If two brain regions show similarities in their BOLD timeseries, they are functionally connected




Types of connectivity

- Functional connectivity
 - Statistical dependency
- Dynamic connectivity
 - Changes in functional connectivity over time
- Effective connectivity
 - Directional influence
- Anatomical (structural) connectivity
 - Presence of a white matter tract






Two broad categories of resting state methods



Voxel-based methods



Node-based methods



Overview of resting state methods

Voxel-based methods

- Seed-based correlation analysis
 - SCA
- Independent component analysis
 - ICA
- Amplitude of low frequency fluctuations
- Regional homogeneity

Node-based methods

- Network modelling analysis
 - FSLNets
- Graph theory analysis
 - Such as degree, hub, path length
- Dynamic causal modelling
- Non-stationary methods



Overview of resting state methods

Voxel-based methods

- Seed-based correlation analysis
 - SCA
- Independent component analysis
 - ICA -> this afternoon
- Amplitude of low frequency fluctuations
- Regional homogeneity

Node-based methods

- Network modelling analysis
 - FSLNets -> tomorrow morning
- Graph theory analysis
 - Such as degree, hub, path length
- Dynamic causal modelling
- Non-stationary methods




Resting state fMRI and ICA

- Introduction to resting state
- Independent Component Analysis
- Single-subject ICA
- Multi-subject ICA
- Dual regression




Resting state functional components


Spatial Maps or Configuration



Time Courses



Network Modelling




- Each component -> a set of brain regions that work in synchrony
- Timecourses can then be used for network modelling



MELODIC: ICA tool in FSL


Melodic

multivariate linear decomposition:




Melodic

multivariate linear decomposition:




Melodic

multivariate linear decomposition:



Melodic


multivariate linear decomposition:



Data is represented as a 2D matrix and decomposed into components



What are components?



- express observed data as linear combination of spatio-temporal processes.
- Components can be liberally recombined to reconstruct the original data.



Some characteristics of MELODIC

- Multivariate voxel-based approach
- Exploratory “model-free” method to find interesting structure in the data
- Gives “spatially independent” components
- “Avoiding overfitting” through automatic model order selection
- “Thresholding” to remove background from main signal



Some characteristics of MELODIC


- Multivariate voxel-based approach
- Exploratory “model-free” method to find interesting structure in the data
- Gives “spatially independent” components
- “Avoiding overfitting” through automatic model order selection
- “Thresholding” to remove background from main signal



ICA vs GLM: Exploratory vs Confirmatory

Melodic

multivariate linear decomposition:




Data is represented as a 2D matrix and decomposed into components



Melodic

multivariate linear decomposition:



Data is represented as a 2D matrix and decomposed into components




Model-based (GLM) analysis

A diagram illustrating the General Linear Model (GLM) equation. It shows a vertical vector of jagged lines representing the measured time-series on the left. This is followed by an equals sign, then a vertical vector of smoother, wavy lines representing the expected response, followed by the Greek letter beta with a subscript 'i' representing the parameter vector, and finally a plus sign and another vertical vector of jagged lines representing the noise term on the right.

- Model each measured time-series as a linear combination of signal and noise
- We know the expected response -> use that to define the design matrix



Resting state = Model-free?



Resting state timeseries are unconstrained -> no design matrix



Task: Confirmatory vs Rest: Exploratory

Confirmatory

- “How well does my model fit to the data?”

Problem \Rightarrow Data \Rightarrow

Model \Rightarrow Analysis

\Rightarrow Results

- results depend on the model

Exploratory

- “Is there anything interesting in the data?”

Problem \Rightarrow Data \Rightarrow

Analysis \Rightarrow Model

\Rightarrow Results

- can give unexpected results



Some characteristics of MELODIC


- Multivariate voxel-based approach
- Exploratory “model-free” method to find interesting structure in the data
- Gives “spatially independent” components
- “Avoiding overfitting” through automatic model order selection
- “Thresholding” to remove background from main signal



Spatially independent components



Spatial ICA for FMRI



- data is decomposed into a set of **spatially independent** maps and a set of time-courses








PCA vs. ICA ?

Simulated
Data

(2 components, slightly
different timecourses)








PCA vs. ICA ?

Simulated
Data

(2 components, slightly
different timecourses)







PCA vs. ICA ?



Simulated Data

(2 components, slightly different timecourses)



PCA



- Timecourses orthogonal
- Spatial maps and timecourses “wrong”







Simulated Data

(2 components, slightly different timecourses)





PCA

- Timecourses orthogonal
- Spatial maps and timecourses “wrong”




ICA

- Timecourses non-co-linear
- Spatial maps and timecourses “right”



PCA vs. ICA


- PCA finds projections of maximum amount of variance in Gaussian data (uses 2nd order statistics only)



Gaussian data

PCA vs. ICA

- PCA finds projections of maximum amount of variance in Gaussian data (uses 2nd order statistics only)
- Independent Component Analysis (ICA) finds projections of maximal independence in non-Gaussian data (using higher-order statistics)




non-Gaussian
data



Non-Gaussianity




sources




mixtures



Non-Gaussianity



non-Gaussian



Gaussian



ICA estimation

- **Random** mixing results in **more** Gaussian-shaped PDFs (Central Limit Theorem)
- ICA turns this around:

if we estimate components with **less** Gaussian distributions this is unlikely to be a random result

➡ measure non-Gaussianity
- can use **neg-entropy** as a measure of non-Gaussianity





ICA estimation

- Need to find an **unmixing matrix** such that: a) maximises independence between components and b) maximises non-Gaussianity of components
- For (a) need a **contrast (objective/cost) function** to drive the unmixing which measures statistical independence and for (b) need an **optimisation technique**:
 - gradient descent & kurtosis or cumulants (**Jade**)
 - gradient descent & maximum entropy (**Infomax**)
 - fixed point iteration & neg-entropy (**FastICA**)



Some characteristics of MELODIC

- Multivariate voxel-based approach
- Exploratory “model-free” method to find interesting structure in the data
- Gives “spatially independent” components
- “Avoiding overfitting” through automatic model order selection
- “Thresholding” to remove background from main signal




Managing overfitting: Automatic model order selection







The 'overfitting' problem

Fitting a noise-free model to noisy observations:

- No control over signal vs. noise (non-interpretable results)
- Statistical significance testing not possible



GLM analysis




standard ICA (unconstrained)



Probabilistic ICA model in MELODIC

statistical “latent variables” model: we observe linear mixtures of hidden sources in the presence of Gaussian noise




Issues:

- Model Order Selection: how many components?




Model Order Selection


‘How many components’?



under-fitting: the amount of explained data variance is insufficient to obtain good estimates of the signals




over-fitting: the inclusion of too many components leads to fragmentation of signal across multiple component maps, reducing the ability to identify the signals of interest



optimal fitting: the amount of explained data variance is sufficient to obtain good estimates of the signals while preventing further splits into spurious components



Model Order Selection



- observed Eigenspectrum of the data covariance matrix
- - - Laplace approximation of the posterior probability of the model order
- . - theoretical Eigenspectrum from Gaussian noise




Some characteristics of MELODIC

- Multivariate voxel-based approach
- Exploratory “model-free” method to find interesting structure in the data
- Gives “spatially independent” components
- “Avoiding overfitting” through automatic model order selection
- “Thresholding” to remove background from main signal




Thresholding

Thresholding




Thresholding

- classical null-hypothesis testing is invalid
- data is assumed to be a linear combination of signals and noise
- the distribution of the estimated spatial maps is a mixture distribution!



Alternative Hypothesis Test



- use Gaussian/Gamma mixture model fitted to the histogram of intensity values (using EM)



Summary of part 1

- Resting state allows us to study the intrinsic organisation of the brain
- It focuses on connectivity and estimates functional components (resting state networks, RSN)
- Each RSN is characterised with a spatial map and a time course
- ICA can be used to characterise these RSNs
- MELODIC is FSL's ICA tool
 - Model-free, spatial independence, non-Gaussianity, model order selection, thresholding via mixture modelling



Resting-State fMRI: ICA and Dual Regression

Rezvan Farahibozorg

FSL Course 2024

19 June, Osaka, Japan



Resting state fMRI and ICA


- Introduction to resting state
- Independent Component Analysis
- **Single-subject ICA**
- Multi-subject ICA
- Dual regression



The goal of single subject ICA is
artefact detection

FMRI inferential path

Experiment



Physiology




Image from Servier Medical Art, released under the CC BY 3.0 license

MR Physics




Image from mos.ru, released under the CC BY 4.0 license

Analysis



Interpretation
of final results






Variability in fMRI

Experiment




suboptimal event timing,
inefficient design, etc.

Physiology




secondary activation, ill-
defined baseline, resting-
fluctuations etc.

MR Physics




MR noise,
field inhomogeneity,
MR artefacts etc.

Analysis



filtering & sampling artefacts, design
misspecification, stats &
thresholding issues etc.

Interpretation
of final results






Artefact detection

- Resting state fMRI data contain a variety of source processes
- Some of these sources are interesting signals and others are artefacts such as motion, cardiac pulsation, respiration
- Artifactual sources typically have unknown spatial and temporal extent and cannot easily be modelled accurately
- Good news: they are spatiotemporally distinct from true signals in the brain
- Therefore, ICA which is an exploratory tool which requires minimal knowledge of spatiotemporal characteristics of components can be used to identify artefacts, a.k.a, **Structured noise**.
- This is done at **single-subject level**.




Manual labelling of ICA components as
signal vs artefact

FSLeyes Melodic Mode




motion





cardiac






susceptibility motion






multiband






signal




effects of scan parameters





manual classification







Removing artefacts




ICA-based denoising






ICA-based denoising





ICA-based denoising






Semi-automatic artefact detection



Semi-automatic classification





semi-automatic classification

- **FIX** (https://fsl.fmrib.ox.ac.uk/fsl/docs/#/resting_state/fix)
 - Classifier with many features
 - Requires manually labelled training data
 - 99% accuracy on high-quality data
- The latest version is reimplemented in Python and installed as part of FSL, unlike previous MATLAB/R version that needed separate installation.



semi-automatic classification

- **FIX** (https://fsl.fmrib.ox.ac.uk/fsl/docs/#/resting_state/fix)
 - Classifier with many features
 - Requires manually labelled training data
 - 99% accuracy on high-quality data
- **ICA-AROMA** (github.com/rhr-pruim/ICA-AROMA)
 - Simple classifier with only 4 features
 - No training data required
 - Mainly designed for motion artefacts



Resting state fMRI and ICA

- Introduction to resting state
- Independent Component Analysis
- Single-subject ICA
- Multi-subject ICA
- Dual regression



The goal of multi subject ICA is to characterise Resting State Networks (RSNs)





Different ICA models

Single-Session ICA

each ICA component comprises:



spatial map & timecourse





Multi-Session or Multi-Subject ICA: Concatenation approach

each ICA component comprises:



spatial map & timecourse

(that can be split up into subject-specific chunks)






Multi-Session or Multi-Subject ICA: Tensor-ICA approach

each ICA component comprises:



spatial map, session-long-timecourse
& subject-strength plot






Different ICA models


Single-Session ICA

each ICA component comprises:
spatial map & timecourse




Multi-Session or Multi-Subject ICA: Concatenation approach

good when:
each subject has **DIFFERENT** timeseries
e.g. resting-state FMRI



Multi-Session or Multi-Subject ICA: Tensor-ICA approach


good when:
each subject has **SAME** timeseries
e.g. activation FMRI





Concatenated ICA


- Concatenate all subjects' data temporally
- Then run ICA
- More appropriate than tensor ICA (for RSNs)





Concatenated ICA

- Data sets must be registered to a common space (anatomical alignment)
- Memory optimisation trick (called MIGP)







Resting state multi-subject ICA

- Why not just run ICA on each subject separately?
- Correspondence problem (eg RSNs across subjects)
- Different splittings sometimes caused by small changes in the data (naughty ICA!)
- Instead - start with a “group-average” ICA
 - But then need to relate group maps back to the individual subjects





Resting state networks

Default Mode Network



Dorsal Attention Network





Resting state fMRI and ICA

- Introduction to resting state
- Independent Component Analysis
- Single-subject ICA
- Multi-subject ICA
- Dual regression




Resting state multi-subject ICA

Group ICA map



Example subject maps derived from dual regression






Dual Regression


Two steps that both involve multiple regression:

1. Extract subject timeseries
2. Extract subject maps





Dual Regression







Dual Regression

Group ICA map



Example subject maps derived from dual regression





Running dual_regression

```
beckmann — bash — bash — 142x23

[islay:~] dual_regression.sh

dual_regression v0.5 (beta)

***NOTE*** ORDER OF COMMAND-LINE ARGUMENTS IS DIFFERENT FROM PREVIOUS VERSION

Usage: dual_regression <group_IC_maps> <des_norm> <design.mat> <design.con> <n_perm> <output_directory> <input1> <input2> <input3> .....
e.g.  dual_regression groupICA.gica/groupmelodic.ica/melodic_IC 1 design.mat design.con 500 grot `cat groupICA.gica/.filelist`

<group_IC_maps_4D>      4D image containing spatial IC maps (melodic_IC) from the whole-group ICA analysis
<des_norm>              0 or 1 (1 is recommended). Whether to variance-normalise the timecourses used as the stage-2 regressors
<design.mat>            Design matrix for final cross-subject modelling with randomise
<design.con>            Design contrasts for final cross-subject modelling with randomise
<n_perm>                Number of permutations for randomise; set to 1 for just raw tstat output, set to 0 to not run randomise at all.
<output_directory>     This directory will be created to hold all output and logfiles
<input1> <input2> ...   List all subjects' preprocessed, standard-space 4D datasets

<design.mat> <design.con> can be replaced with just
-1                      for group-mean (one-group t-test) modelling.
If you need to add other randomise option then just edit the line after "EDIT HERE" below


[islay:~] █
```

- FSL command line tool, combining:
 - DR to create subject-wise estimates (stage 1 + stage 2)
 - Group comparison using randomise (stage 3)






Group comparison

- Collect maps and perform voxel-wise test (e.g. randomisation test on GLM)



- Can now do voxelwise testing across subjects, separately for each original group ICA map
- Can choose to look at strength-and-shape differences

Group analysis on maps





- can use the Glm tool (Glm_gui on mac) to create GLM design and contrast matrices



Dual regression outputs

- `dr_stage1_subject[#SUB].txt` - the timeseries outputs of stage 1 of the dual-regression.
- `dr_stage2_subject[#SUB].nii.gz` - the spatial maps outputs of stage 2 of the dual-regression.
- `dr_stage2_ic[#ICA].nii.gz` - the re-organised parameter estimate images
- `dr_stage3_ic[#ICA]_tstat[#CON].nii.gz` - the output from randomise




(corrected for multiple comparisons across voxels
but not across #components!!)



Group template maps

- Generate from the data using ICA
 - use all data to get unbiased templates
 - use independent control group
 - will not model signals and artefacts
- use existing template



template RSNs

<https://www.fmrib.ox.ac.uk/datasets/royalsoc8/>



Summary of part 2

- MELODIC ICA is applied to resting state fMRI for two main purposes: artefact detection, resting state network characterisation
- Artefact detection is done using single subject ICA
 - Artefact components can be classified manually or using FIX/AROMA
 - We then use GLM to regress out the artefacts and clean up the fMRI
- RSN estimation is done using group-level ICA
 - We use concatenation technique
 - This will give consensus RSNs which we then project onto subject fMRI using dual regression
 - Subject specific RSNs obtained from dual regression can be used for statistical comparisons (and input to FSLNets -> we will learn about this tomorrow)



Resting state fMRI and ICA

Available from:

- [Oxford University Press](#)
- [Amazon](#)

